



## ADAPTATIONS FOR TOILETING

Reference: Awareness-Cum-Training Packages In Disability, Training In Daily Living Skills And Managing Difficult Behaviours. Unit-5, IGNOU School Of Continuing Education




Cushion or board to cover the hole in the wheelchair


Adapting the wheelchair for toileting



A hole carved in the seat of the chair, which can be placed over the w.c. During the rest of the day, the hole can be covered with a wooden plank, with a cushion over it, and the chair can be used for sitting.



Various types of potties that provide varying levels of back support



A wooden frame fixed around the potty provides back and side support to the child



Bars inside the toilet help the child to maintain balance while squatting